





## Academia Sinica Protein Clinic Report


Date	2020.01.06
Case number	20200106_test_1
Applicant	Hsien-Ya Lin
Principle investigator	Chun-Hung Hans Lin
Affiliation	Academia Sinica
Institute/department	Institute of Biological Chemistry
ASPC project manager	<input checked="" type="checkbox"/> Hsien-Ya Lin <input type="checkbox"/> Yun-Ru Chen
Number of sample(s)	1
<b>Sample 1</b>	
Sample name	<i>Bifidobacterium dentium</i> $\beta$ -glucuronidase
Number of amino acids	692
Molecular weight (Da)	77065
Theoretical pI	4.99
Extinction coefficient @280 nm ( $M^{-1} cm^{-1}$ )	119180
Buffer condition	150mM NaCl, 50mM Tris-HCl (pH 8.0), 0.02% NaN <sub>3</sub>
Package option	<input checked="" type="checkbox"/> Deluxe <input type="checkbox"/> Basic <input type="checkbox"/> SEC-MALS only
Service items	<input checked="" type="checkbox"/> CF-01 SDS-PAGE
	<input checked="" type="checkbox"/> CF-02 UV-Vis scan & 260/280
	<input checked="" type="checkbox"/> CF-03 Tycho NT. 6
	<input checked="" type="checkbox"/> CF-04 Far-UV CD spectrum
	<input checked="" type="checkbox"/> CF-05 DSF T <sub>m</sub> screen
	<input checked="" type="checkbox"/> CF-06 AUC (SV mode)
	<input type="checkbox"/> CF-07 SEC
	<input checked="" type="checkbox"/> CF-08 SEC-MALS
Supporting information	<input checked="" type="checkbox"/> Tycho NT. 6 raw data
	<input checked="" type="checkbox"/> Far-UV CD raw data
	<input checked="" type="checkbox"/> DSF screen condition details
	<input checked="" type="checkbox"/> AUC raw data
	<input checked="" type="checkbox"/> SEC-MALS Wyatt report

*Disclaimer: The content is for research purpose only. Do not use this report for commercial purpose. ASPC only performs the analyses. ASPC is not responsible for the sample preparation.*


The content is for research purpose only. Do not use this report for commercial purpose.

CF-01 SDS-PAGE		CF-02 UV-Vis scan and 260/280 ratio	
SDS-PAGE	4-20% gradient	Instrument	Beckman DU <sup>®</sup> 730
Protein marker	Thermo Scientific PageRuler™ Prestained Protein Ladder, 10 to 180 kDa	Sample concentration (μM)	3
Protein sample Protein:4x loading dye=3.1	5 μL	OD <sub>260</sub> /OD <sub>280</sub> *Ration < 0.7: no nucleic acid contamination	0.57
			
<b>CF-03 Tycho NT.6</b>			
The system measures the fluorescence of intrinsic tryptophan and tyrosine residues detected at both 350 nm and 330 nm as a 30°C/min. temperature ramp is applied from 35–95°C.			
<p><b>Profiles Ratio 350 nm / 330 nm</b></p> 		<p><b>Profiles First derivative (ratio)</b></p> 	
Inflection point (°C)			<b>67.2</b>



The content is for research purpose only. Do not use this report for commercial purpose.



The content is for research purpose only. Do not use this report for commercial purpose.



The content is for research purpose only. Do not use this report for commercial purpose.

<b>CF-06 AUC SV mode</b>	
Data collection parameters:	
Sample buffer	150mM NaCl, 50mM Tris pH8.0, 0.02% NaN <sub>3</sub>
Detection mode	UV absorption at 280nm
Rotation speed (rpm)	40 000
Duration (hour)	5
Sample temperature (°C)	20
Buffer density (g/ml)*	1.0367
Buffer viscosity (η)*	0.0105373
* calculated by Sednterp based on the buffer composition	
Raw data fitting	
 <p>The figure displays raw data fitting for sedimentation velocity. The top panel shows signal (AU) on the y-axis (0.0 to 0.5) versus radius (cm) on the x-axis (6.2 to 7.0). Multiple sigmoidal curves are shown, representing different time points. The bottom panel shows residuals on the y-axis (-0.03 to 0.02) versus radius (cm) on the x-axis (6.2 to 7.0), indicating the fit quality.</p>	
SV distribution	
 <p>The sedimentation velocity distribution plot shows c(s) (AU/S) on the y-axis (0.0 to 0.6) versus sedimentation coefficient (S) on the x-axis (0 to 20). A major peak is observed at approximately 11.5 S, with smaller peaks at approximately 4.5 S and 19.5 S.</p>	
Molecular weight estimate (kDa)	84.5 (monomer), 150 (dimer) and 391 (tetramer)

The content is for research purpose only. Do not use this report for commercial purpose.

<b>CF-07 SEC-MALS</b>	
Data collection parameters:	
Sample buffer	150mM NaCl, 50mM Tris pH8.0, 0.02% NaN <sub>3</sub>
SEC column	GE 10/300 Superdex 200 increase
Flow rate	0.5 ml/min
Sample temperature (°C)	ambient
Control sample (BSA)	1.0 mg/ml
Calculated MW of BSA (kDa)	131.9 ± 0.82 (dimer)
Theoretical value 66.5 kDa	68.09 ± 0.32 (monomer)
<div style="text-align: center;">Results Fitting</div>	
Sample 1	1.0 mg/ml
Calculated MW (kDa)	(peak 1) 308.5 ± 0.10 (tetramer) (peak 2) 167.6 ± 1.43 (dimer) (peak 3) 82.13 ± 10.84 (monomer)
Polydispersity	(peak 1) Mw/Mn: 1.000 (± 0.529%), Mz/Mn: 1.000 (± 0.917%) (peak 2) Mw/Mn: 1.009 (± 15.108%), Mz/Mn: 1.019 (± 27.180%) (peak 3) Mw/Mn: 1.865 (± 92.616%), Mz/Mn: 3.163 (± 175.173%)

The content is for research purpose only. Do not use this report for commercial purpose.

