





## Academia Sinica Protein Clinic Report


Date	2020.01.06
Case number	20200106_test_2
Applicant	Yun-Ru Chen
Principle investigator	Shang-Te Danny Hsu
Affiliation	Academia Sinica
Institute/department	Institute of Biological Chemistry
ASPC project manager	<input type="checkbox"/> Hsien-Ya Lin <input checked="" type="checkbox"/> Yun-Ru Chen
Number of sample(s)	1
<b>Sample 1</b>	
Sample name	BAP1
Number of amino acids	238
Molecular weight (Da)	26949.72
Theoretical pI	5.63
Extinction coefficient @280 nm ( $M^{-1} cm^{-1}$ )	30940
Buffer condition	PBS pH7.4, 10% glycerol, 0.02% NaN <sub>3</sub>
Package option	<input checked="" type="checkbox"/> Deluxe <input type="checkbox"/> Basic <input type="checkbox"/> SEC-MALS only
Service items	<input checked="" type="checkbox"/> CF-01 SDS-PAGE
	<input checked="" type="checkbox"/> CF-02 UV-Vis scan & 260/280
	<input checked="" type="checkbox"/> CF-03 Tycho NT. 6
	<input checked="" type="checkbox"/> CF-04 Far-UV CD spectrum
	<input checked="" type="checkbox"/> CF-05 DSF T <sub>m</sub> screen
	<input checked="" type="checkbox"/> CF-06 AUC (SV mode)
	<input type="checkbox"/> CF-07 SEC
	<input checked="" type="checkbox"/> CF-08 SEC-MALS
Supporting information	<input checked="" type="checkbox"/> Tycho NT. 6 raw data
	<input checked="" type="checkbox"/> Far-UV CD raw data
	<input checked="" type="checkbox"/> DSF screen condition details
	<input checked="" type="checkbox"/> AUC raw data
	<input checked="" type="checkbox"/> SEC-MALS Wyatt report

*Disclaimer: The content is for research purpose only. Do not use this report for commercial purpose. ASPC only performs the analyses. ASPC is not responsible for the sample preparation.*


The content is for research purpose only. Do not use this report for commercial purpose.

CF-01 SDS-PAGE		CF-02 UV-Vis scan and 260/280 ratio	
SDS-PAGE	15% gradient	Instrument	Beckman DU® 730
Protein marker	Thermo Scientific PageRuler™ Prestained Protein Ladder, 10 to 180 kDa	Sample concentration (μM)	24.14
Protein sample Protein:4x loading dye=3.1	5 μL	OD <sub>260</sub> /OD <sub>280</sub> *Ration < 0.7: no nucleic acid contamination	0.66
			
<b>CF-03 Tycho NT.6</b>			
The system measures the fluorescence of intrinsic tryptophan and tyrosine residues detected at both 350 nm and 330 nm as a 30°C/min. temperature ramp is applied from 35–95°C.			
			
		Inflection point (°C)	<b>52.4</b>



The content is for research purpose only. Do not use this report for commercial purpose.




The content is for research purpose only. Do not use this report for commercial purpose.




The content is for research purpose only. Do not use this report for commercial purpose.

<b>CF-06 AUC SV mode</b>	
Data collection parameters:	
Sample buffer	PBS pH7.4, 10% glycerol, 0.02% NaN <sub>3</sub>
Detection mode	UV absorption at 280nm
Rotation speed (rpm)	40 000
Duration (hour)	5
Sample temperature (°C)	20
Buffer density (g/ml)*	1.034
Buffer viscosity (η)*	0.0138995
* calculated by Sednterp based on the buffer composition	
Raw data fitting	
 <p>The plot displays multiple sedimentation curves for different time points, showing signal intensity (AU) as a function of radius (cm). The curves are color-coded and show a characteristic sigmoidal shape. Below the main plot is a residuals plot showing the difference between the experimental data and the fitted curves, which fluctuates around zero.</p>	
SV distribution	
 <p>The sedimentation velocity (SV) distribution plot shows a single, sharp peak at a sedimentation coefficient of approximately 1.5 S. The y-axis represents the concentration c(s) in AU/S, ranging from 0.0 to 3.5. The x-axis represents the sedimentation coefficient in S, ranging from 0 to 10.</p>	
Molecular weight estimate (kDa)	26.6 (monomer)

The content is for research purpose only. Do not use this report for commercial purpose.



The content is for research purpose only. Do not use this report for commercial purpose.



CONFIDENTIAL